Instructions for use: ELITEsmile Light Whitening 32% Hydrogen Peroxide (mixed) - medical device

Tooth Whitening System Information for the dentist – Medical tooth whitening

Medical Tooth Whitening System for use in the dental clinic for strongly discoloured vital and non-vital teeth, only with medical indications. In accordance with the Medical Device Directive 93/42/EEC. ELITEsmile Light Whitening 32% (mixed) complies with standard EN ISO 28399.

ELITEsmile LIGHT WHITENING Gel is available in a mixed concentration of 32% Hydrogen Peroxide recommended for the chair side use with whitening lamps

Storage and shelf life

ELITEsmile Light Whitening 32% (MIXED) gel should be stored in a cool place or refrigerator (4 – 12 °C).

The expiry date on the package or syringe refers only to products that have been stored correctly. Do not freeze the gel and protect it from heat and sunlight.

Keep out of the reach of children. Do not use after expiry date.

IMPORTANT NOTE: Please carefully read and understand the instructions before proceeding with tooth the whitening treatment.

Description

ELITEsmile Light Whitening 32% (MIXED) is a 2-component material for whitening discolored teeth. It is intended only for professional medical whitening in the dental office. By mixing the two components of the gel, the bleaching agent is activated. The duration of the treatment is determined individually and depends on the sensitivity of the teeth, to the desired degree of whitening and the nature and extent of tooth discolorations.

The colored gel can be applied directly from the mixing tip to the teeth - its viscous consistency prevents accidental dripping or flow.

Unused gel can be kept until the end of the shelf life in the closed syringe. It can be used, along with a new mixing tip, for another treatment.

Ingredients

Aqua, Hydrogen Peroxide, Silica, Glycerol, Organic amines, Polyglycols, coloring agent

Indications (fields of application)

ELITEsmile Light Whitening 32% (MIXED) is used for fast, medically indicated whitening of strongly discolored, vital and non-vital teeth, which cannot be converted back by a professional dental cleaning into the natural or normal tooth color. Medically tooth whitening is indicated when discolorations lead to a deterioration of (health- related) quality of life of the patient due to their abnormal deviating tooth color. Medical tooth whitening can also be used as an alternative to (significantly more invasive) conservative treatment of the affected teeth with veneers or crowns.

Tooth colors, as they normally appear in their brightness, translucency and color in the population, are considered natural or normal. For diagnosis and indication e.g. the A1-D4 Vita classical shade guide can be used (for the most frequently occurring tooth colors).

Examples of the causes of such (intrinsic) discoloration are: tetracycline/minocycline discolorations, moderate to severe dental fluorosis, non-vital teeth (pulp necrosis) – e.g. by trauma, genetic or age-related (degenerative) discoloration, porphyria, or jaundice, heavy deposits of foods with penetration into the tooth substance (intrinsic discoloration by extrinsic dyes).

It can be used to whiten single or multiple teeth or parts of teeth. It can also be used for whitening of single endodontically treated non-vital teeth including intra-coronal bleaching.

Application

Prior to treatment

Before the bleaching treatment a diagnostic indication for "medical indicated tooth whitening" has to be defined by the dentist. After examining the patient, the current tooth color should be documented. Professional tooth cleaning is indicated prior (ideally 1-2 weeks) to the whitening treatment.

Because restorative materials will not whiten, it is recommended that teeth be whitened prior to placing aesthetic restorations. Any dental work (crowns, fillings etc.) should be done 1 – 2 weeks after the treatment in order to match the color of the newly whitened teeth. See also Interactions.

1) Start the whitening procedure by applying a thin layer of lip protection on the patient's lips. The teeth should be kept free, by using a lip retractor (not provided in the kit, can be ordered separately), in the patient's mouth. If necessary, use cotton rolls.

2.) Application of the GINGIVA PROTECTOR

GINGIVA PROTECTOR is an easy to apply light-curing resin, which allows optimal covering of soft tissue and neighboring teeth during whitening. It is also suitable to cover exposed tooth necks and filling materials during the treatment.

Attach the dispensing tip to the syringe. Dry the gingiva carefully to ensure that the GINGIVA PROTECTOR will hold and seal during the whole treatment.

Apply the GINGIVA PROTECTOR direct on the gingiva and cover 3-4 mm of the gingiva and 0,5-1 mm of the teeth. Work in quadrants: light cure the gingiva protector after each half jaw finished in an arching movement for about 25 seconds. When finished, use dental mirror to check if the seal is tight and no reddish gingiva is visible. Close open spots again with Gingiva Protector and light cure.

Instead of the GINGIVA PROTECTOR a rubber dam can also be used.

3.) Whitening procedure

The gel should be used at room temperature. Attach the mixing tip to the top of the dual-syringe of the ELITEsmile Light Whitening 32% (MIXED) gel.
Verify the flow of the gel on gauze or mixing pad, prior to the intraoral application.

Apply the gel directly to the surface of the teeth. The gel should be applied about 1 – 2 mm thick (avoid contact with other parts of the body and clothes).

The patient should be positioned in a slightly declined position (45°-degree angle) during the whitening process to be able to swallow saliva comfortable.

Place a whitening lamp in front of the patient's mouth during the gel application and turn on the lamp. Eye protection with safety goggles for patient/operator are mandatory.

Allow the gel to work about 15-20 minutes. Reduce the application time and cycles if the patient experiences sensitivity.

Remove the gel from the teeth using suction – do not rinse – then apply a fresh layer of gel to the teeth. This process can be repeated three times.

(Maximum treatment time is 60 minutes - 3x 15-20min)

After the last treatment remove the gel first using suction and then by rinsing it off by water.

Now the GINGIVA PROTECTOR and lip holder can be removed. Finally, the patient should rinse their mouth with warm water.

The tooth color after the treatment should be documented for your records.

Finish off the treatment and lower the risk of unpleasant sensitivity for the patient by using After Whitening Mousse – apply to the teeth and leave on for about 10 minutes. Recommend WHITESMILE WHITENING TOOTHPASTE with 30% Xylitol to your patient to maintain the results and to protect the teeth.

Internal bleaching of endodontically treated, non-vital teeth

Wait at least one week following the endodontic procedure before performing an internal tooth whitening and ensure that the tooth is endodontically sound.

1. Check the X-rays, to ensure that an adequate apical endodontic sealer is in place and to check where the alveolar bone margin runs.
2. Remove the filling material from the access hole until the depth you defined on basis of the X-ray image. Continue to remove the filling material and / or temporary endodontic materials in order to allow the tooth whitening agent to get in contact with the internal structure of the tooth for best results.
3. Seal the root with conventional glassionomer cement. Apply a 1-2 mm thick layer to ensure a secure seal.
4. The internal whitening is performed according to the instructions for vital teeth. Cover the gingiva around the tooth with Gingiva Protector. Apply ELITEsmile Light Whitening 32% (MIXED) gel into the cavity and make sure the tooth structure is covered with 1-2mm gel. Apply the gel also on the outside surface of the tooth. Wait 15-20 minutes. Suck of the gel. Procedure might be repeated for 1 or 2 times. Do not use light activation inside the cavity.

After the last application suck of gel and rinse off remaining gel from cavity.

5. Use non-eugenol temporary restorative material to seal the opening.
6. Before starting to work on adjacent teeth, it is recommended to wait at least two weeks after the whitening treatment - in order to stabilize the tooth color.

Side effects

Sensitivity and irritation of the teeth with more or less pronounced pain during or after treatment can occur in individual cases. As a rule, these disappear within a short time. In case of excess sensitivity of the teeth during treatment, the exposure time of the bleaching agent must be shortened, or treatment discontinued. Nonetheless, to obtain an acceptable bleaching result, fractionated whitening should be conducted over several sessions in such cases.

Sometimes after the whitening treatments white spots appear on the teeth. These spots are mostly from fluorosis that has not been visible before. After a period of a few days these white spots will usually disappear.

Some studies refer to a temporary change in the organic dental matrix as well as a reduction in the microhardness of enamel and dentine, however without any pathological relevance.

In the case of endodontic bleaching, an inadequately performed root canal filling which does not reach or exceed the enamel cement border, may lead to damage of the periodontium (external root resorption) or even loss of the tooth.

It is essential for gingival protection to be extended over the tooth neck (approx. 1 mm). Inadequately applied gingival protection can lead to irritation of the gingiva through contact with peroxide. In this case the bleaching therapy must be discontinued, and the tooth prepared anew after thorough rinsing and drying.

Contains small amounts of organic amines. These can cause allergic reactions when in contact with skin or mucosa.

Sensitivity of the teeth can be reduced with a subsequent application of AFTER WHITENING MOUSSE.

Contraindications

ELITEsmile Light Whitening 32% (MIXED) gel should not be used in the following cases: stains which can be removed by suitable polishing agents and rotating brush during a professional tooth cleaning, Pregnancy and lactation due to lack of studies, known allergy to any of the ingredients, exposed dentin, faulty edges of fillings, unmet primary and secondary caries, very heavy smokers and alcoholics, persons under the age of 18, under anesthesia.

Interactions

Hydrogen peroxide weakens the polymerization of polymers. Therefore, composite fillings or bracket bonds should not be performed earlier than 2 weeks after a whitening treatment, due to lack of adhesion.

Special notes

For at least 4 hours after whitening session, the patient should avoid consuming tea, coffee, red wine, cigarettes etc. as this will reduce effectiveness of treatment.

The bleaching effect can intensify for a few days after completing treatment. The result of bleaching treatment cannot be guaranteed as tooth discoloration may be due to a variety of causes. The bleaching effect can change a few shades in a period of a week after completing treatment.

Any existing composite fillings are not bleached. The earliest these can be renewed, is two weeks after completing bleaching treatment.

Orthodontic work (bonding brackets) should either be performed prior to bleaching therapy or no earlier than 2 weeks after (impaired adhesion to teeth).

Warnings

Because of the high concentration of Hydrogen Peroxide ELITEsmile Light Whitening 32% (MIXED) gel is highly reactive. To avoid contact with skin or eyes we recommend the use of gloves and safety goggles. If stored improperly, the syringes may become pressurized and lead to uncontrolled spillage of gel when opening.

For safety reasons, the patient should also wear safety goggles. If the gel does come into contact with skin, remove it and rinse immediately with water. White spots may appear in the affected area and the patient may experience a burning sensation. The symptoms should disappear after approx. 30 – 60 minutes. If the gel comes into contact with the eyes, rinse immediately and consult a doctor.

Notes on disposal

Full syringes must be sent for special disposal in accordance with regulations on special waste.

This product was developed for use in dentistry and must be used in accordance with the instructions for use. The manufacturer cannot be made liable for any alternative use. In addition, the user is obliged to check the material at his own responsibility for its suitability and application for the intended purposes, especially if these are not listed in the instructions for use.

Symbols

- Batch description 
- Expiry date 
- Observe instructions for use 
- Recommended storage temperature 
- Manufacturer 
- Harmful if swallowed or if inhaled 

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