with 1–2 mm gel.

- In addition, apply whitening gel to the outer surface of the tooth. Leave on for 15–20 minutes, then suction off. Repeat the process 3–5 times.

- After the last application, suction off the whitening gel and rinse the cavity thoroughly.
- Seal the cavity with an eugenol-free temporary restoration.
- Final closure of the canal after approx. 2 weeks, see section Interactions.

Application notes

The whitening result may change by a few shades in the days following the treatment. The result of tooth whitening cannot be guaranteed, as different patients have different causes. Restorative materials such as composite fillings are not whitened. Lack of or inadequate gum protection can lead to irritation and sensitization. One mixing tip can be used for all cycles within a session. The product is intended for single patient use only and should be used as intended. After treatment, existing white spots, e.g. in the case of fluorosis, may become more visible; this is usually a temporary effect. The duration of treatment depends on the sensitivity of the teeth, the desired degree of whitening and the type and extent of discoloration.

Warnings

Serious incidents in connection with this medical device must be reported to BeViA Dental and the competent authority. The manufacturer is not liable for misuse. Please always observe the current safety data sheet.

Tooth whitening gel

The colorless component of the double-chamber syringe is strongly acidic; the second component is strongly alkaline. When unmixed, both are highly corrosive. Always use the correct mixing ratio.

In case of contact with eyes: Rinse carefully with plenty of water for several minutes. Remove contact lenses if possible. Continue rinsing. Consult a doctor.

If swallowed: Rinse out mouth. DO NOT induce vomiting.

In case of contact with skin: Wash carefully with plenty of soap and water. Moisturize with a cream.

Causes serious eye damage. Wear protective gloves / protective clothing / eye protection

Gingiva Protector

In case of skin irritation or rash: Seek medical advice/attention! Allergic potential due to acrylates.

May cause allergic skin reactions. Harmful to aquatic organisms, with long lasting effects.

Storage and shelf life

Store the gel lightly closed at +4–12 °C.

Seal packs as indicated. Protect from heat and sunlight and do not freeze. Avoid the use of dental fluorides. Do not use after the expiry date. Improper storage may result in pressure build-up in the syringes and the components may leak out when opened.

Disposal instructions

The gel must be disposed of in accordance with the applicable requirements for hazardous waste. If the syringes have not been completely squeezed out, close them again with the sealing cap and dispose of them separately. Never re-use the double-chamber syringe or transfer into the waste container and disposal. The resulting heat can ignite other waste components, e.g. paper, and cause a fire. Gingiva Protector must not be disposed of in wastewater.

English

Instructions for use for BMS white 38% HP (30% mixed)

The gel is a 2-component gel for whitening pathologically discolored teeth for use in adult patients (≥ 18 years). It is only intended for professional, medical whitening in dental practices. The product is activated by automatically mixing the two components. The product complies with the EN ISO 23899 standard.

The gel is chemically activated.

Intended purpose

The gel is used for the fast, medicaly indicated whitening of pathologically and severely discolored, vital and non-vital teeth in cases where the natural or normal tooth shade cannot be restored by professional tooth cleaning.

The Gingiva Protector is a light-curing resin designed to protect the gums from chemical burns caused by high-concentration hydrogen peroxide during tooth whitening treatment.

Indications

The gel is used for the fast, medicaly indicated whitening of pathologically and severely discolored, vital and non-vital teeth in cases where the natural or normal tooth shade cannot be restored by professional tooth cleaning. Discolorations are considered medically indicated if they lead to a deterioration in the patient's (oral health-related) quality of life as the tooth shade varies from the norm. If necessary, whitening can also be used as an alternative to covering the affected teeth with veneers or crowns.

Examples of the causes of such (intrinsic) discoloration are tetracycline/metalloids discoloration, moderate to severe dental fluorosis, non-vital teeth (pulp necrosis) – e.g. due to trauma, genetic or age-related (degenerative) discoloration, porphyria or iatrogenic/aural and heavy deposits of food with penetration into the tooth substance by extrinsic pigments.

Tooth shades considered natural or normal are those that are commonly found in the population in terms of brightness, translucency and color. The Vita A1-D4 classical shade guide (for the most frequently occurring natural colors) is used as aid for establishing a diagnosis and indication.

Application is carried out on single or several vital or non-vital teeth. It is also suitable for whitening individual, endodontically treated non-vital teeth, including intracrown whitening.

Target patient group

The product is intended for adult patients (≥ 18 years) with pathological tooth discoloration. It is intended for single use and must not be reprocessed. The duration of treatment depends on the sensitivity of the teeth, the desired degree of whitening and the type and extent of discoloration.

Contraindications

Do not use in patients under 16 years of age, during pregnancy or breast-feeding and in anesthetized patients. Discoloration that can be removed during professional tooth whitening. Application is also excluded in the case of inflamed gingiva, exposed dentin, non-vital teeth, defective fillings or untreated primary or secondary caries. Tooth whitening should be avoided in severely traumatized, tooth-less teeth, as there is an increased risk of undetected crack formation and resulting cervical root caries.

Side effects

Sensitivity or pain may occur during or after the treatment. This usually subsides quickly in the event of a mild irritation. In severe cases, a short-term or discontinue treatment and continue in several sessions if necessary. During endodontic whitening, inadequate canal filling can lead to damage to the periradicular and even tooth loss. Organic amines contained in the product may cause allergic reactions.

The gum protection can cause allergic reactions due to acrylic resins.

Interactions

Hydrogen peroxide can inhibit radical polymerization and thus impair the adhesion of composite fillings or bracket bonding. Such treatments should be carried out no earlier than two weeks after tooth whitening.

Content

Gel: Aqua, Hydrogen Peroxide, Silica, Glycerol, organic amines, polyglycols, coloring agent.

Gingiva Protector: Oligomeric Urethane Acrylate, Medium Chain Triglycerides, Polyacrylic Acid, Difunctional Acrylates, Silica, Zinc Oxide, photoinitiators, pigments.

Application

Tooth whitening should always be carried out before any visible restoration or restorations are placed. Existing restorations should be covered in shade after whitening. For single-tooth whitening, Gingiva Protector can be used to cover teeth that are not to be treated as well as prosthetic restorations (e.g., crowns). During treatment, the patient should be placed in a slightly inclined position (approx. 45°) to facilitate swallowing.

Pre-treatment

- Diagnosis by the dentist before starting treatment.
- Professional tooth cleaning 1–2 weeks before whitening.
- Determine and document the initial tooth shade.
- Protect lips with lip balm, teeth or similar must be kept exposed with the aid of a lip retainer, lip retractor or similar instrument, the patient should wear protective goggles.

- Apply Gingiva Protector: Dry the gingiva with air, attach the applying material in firm, regular quadrants 3–5 mm wide and 2 mm high to the gums and approx. 1 mm to the neck of the tooth, carefully seal the interdental spaces. Cure with a polymerization lamp in an arch-wrapped mode for approx. 25 s (smaller areas 2 x 10 s)

- Check for leaks and seal open areas.

When whitening a single tooth, the gum protection can also be applied to the neighboring teeth so that they are not whitened.

Tooth whitening

The product should be used at room temperature. The mixing tip is placed on the double syringe, a small amount of gel is squeezed out and only homogeneous red gel is used. The patient and practitioner wear protective goggles.

- The whitening gel should be applied in a layer of at least 1 mm on the labial surface of the teeth to be treated.
- Suction off the whitening gel at the end of the cycle (no rinse)
- Further treatment cycles: Repeat steps 1 and 2. The maximum treatment time is 60 min (3 × 20 min), 45 min (3 × 15 min) is recommended. The duration depends on the discoloration, desired degree of whitening and the type and extent of discoloration.

Components

Gel: Aqua, Hydrogen Peroxide, Silica, Glycerol, organic amines, polyglycols, coloring agent.

Gingiva Protector: Oligomeric Urethane Acrylate, Medium Chain Triglycerides, Polyacrylic Acid, Difunctional Acrylates, Silica, Zinc Oxide, photoinitiators, pigments.

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